



SLIT LAMP  
**SL-1E**






#### NOMENCLATURE

- |                                |   |
|--------------------------------|---|
| 1. 10X Eyepieces               | 7. Slit Rotation Control Ring             |
| 2. Chin-Rest Elevation Control | 8. Aperture and Slit Length Control Wheel |
| 3. Slit Width Control Knob     | 9. Filter Selection Lever                 |
| 4. Fixation Target (Option)    | 10. Lamphouse Cover                       |
| 5. Forehead Rest               | 11. Omni-Directional Joystick             |
| 6. Magnification Changer Lever | 12. Cross-Slide Base                      |

#### 14mm Slit Illumination

For maximum flexibility, the slit width is continuously adjustable from 0 to 14mm. The slit length is controlled by seven apertures—1mm, 2mm, 3mm, 5mm, 8mm, 10mm, and 14mm. The 14mm height facilitates the evaluation of today's larger diameter soft contact lenses.



#### Omni-directional Joystick

Fast, effortless positioning is assured by TOPCON's omni-directional joystick. The one-hand control permits precise centering, height adjustment, and focusing.



# atility and Performance You Can Depend On **SL-1E** SLIT LAMP

## Variable Magnification Capability

Comfortable stereoscopic viewing is provided by the converging biomicroscope system. The conveniently placed lever smoothly switches the objective lenses to give 10X and 16X magnification with the 10X eyepieces. The optional 16X eyepieces offer 16X and 25X magnification.



## Built-in Filters

Versatility in examination is enhanced by the three filters which can be easily inserted into the light path. A blue filter for tonometry and fluorescein excitation, a green filter for vascular contrast, and a heat-absorbing filter for patient comfort are incorporated.

## Halogen Illumination

Bright, clean illumination is provided by the halogen lamp. Visualization is improved due to the properties of halogen illumination—the blue characteristic supplies a more natural view, while a darkening of the bulb does not occur through use.



## OPTIONAL ACCESSORIES

### 16X Eyepieces

Allows increased magnification for expanded examinations.

### Hruby Lens

For observation of the fundus or the posterior segment of the vitreous body.



### Fixation Target

The illuminated fixation target can be placed on either side of the chin-rest to assist in eliciting patient cooperation.



### 10X Measuring Eyepiece

When used in place of the standard eyepieces, its millimeter scale permits accurate linear measurements while the protractor scale in 5° increments aids in the fitting of toric contact lenses.

### Applanation Tonometer

The Haag-Streit Model R-900 Tonometer is recommended.

## SPECIFICATIONS

### MICROSCOPE

Type	Binocular Stereoscopic microscope with erect image.		
Magnification changer	Revolver type objective changing system provides two magnification changes		
Objective	1X and 1.6X		
Eyepieces	10X and 16X (Option)		
Magnification & field of view			
Objective × Eyepieces = Magnification	Field of view		
1X 10X 10X	18mm diameter		
1.6X 10X 16X	14.5mm diameter		
1X 16X 16X	11.25mm diameter		
1.6X 16X 25.6X	9mm diameter		
Pupillary adjustments	55mm to 75mm		
Diopter adjustments	- 5D ~ + 3D		

### SLIT ILLUMINATION

Slit projection	1.23X
Slit width	Continuously variable from 0 ~ 14mm
Slit length	14mm, 10mm, 8mm, 5mm, 3mm, 2mm, 1mm
Aperture diameters	14mm, 10mm, 8mm, 5mm, 3mm, 2mm, 1mm
Slit angle	Continuously rotatable from vertical to horizontal
Filters	Blue and Red-free
Lamp	6V 20W halogen

### BASE

Longitudinal movement	80mm
Lateral movement	100mm
Fine cross-slide adjustments	12mm
Vertical movement	30mm

### CHIN-REST

Vertical movements	80mm
--------------------	------

### POWER

Input (primary)	AC 100V, 120V, 220V and 240V; adjustable with built-in voltage selector, 50/60Hz
Output (secondary)	AC 5.5V and 6.3V
Power consumption	30VA

### DIMENSION & WEIGHTS

Unit model	440mm × 350mm 18kgs.
Table model	560mm × 350mm 18kgs.

### OPTIONAL ACCESSORIES

16X Eyepiece	
Fixation target	
Hruby lens	For examinations of the vitreous body and fundus.
10X Measuring Eyepiece	
Accessory drawer	

Subject to change in design and/or specifications without advance notice.

**IMPORTANT**

In order to obtain the best results with this instrument,  
please be sure to review all user instructions prior to operation.



Certificate No.95/4370  
TOPCON CORPORATION  
Ophthalmic & Medical  
Instruments Division (Tokyo)



75-1 Hasunuma-cho, Itabashi-ku, Tokyo 174, Japan. Phone: 3-3558-2520 Fax: 3-3960-4214

**TOPCON AMERICA CORPORATION**

CORPORATE OFFICE  
37 West Century Road, Paramus, NJ 07652, U.S.A.  
Phone: 201-261-9450 Fax: 201-387-2710 WWW: Topcon.com

**TOPCON OMNI SYSTEMS, INC.**

Valley Forge Business Center, 2430 Blvd. of the Generals,  
Norristown, PA 19403, U.S.A.  
Phone: 610-830-9200 Fax: 610-630-6428

**TOPCON EUROPE B.V.**

Esse Basn 11, 2808 LJ Capelle afd Ussel, THE NETHERLANDS  
Phone: 010-4585077 Fax: 010-4585045

**TOPCON S.A.R.L.**

HEAD OFFICE  
104/106, Rue Pilvay, 92300, Levallois-Perret, FRANCE  
Phone: 01-41069494 Fax: 01-47390251

**LYON OFFICE**

138, Avenue du 8 Mai 1945, 69100 Villeurbanne, FRANCE  
Phone: 0478658237 Fax: 0478681902

**TOPCON DEUTSCHLAND GmbH**

Hallesstr. 7, 47677 Willich, GERMANY  
Phone: 02154-9290 Fax: 02154-929111 Telex: 8531981 TOPC D

**TOPCON ESPAÑA S.A.**

HEAD OFFICE  
Frederic Mompou 5, 08960 Sant Just Desvern Barcelona, SPAIN  
Phone: 03-4734057 Fax: 03-4733932

**MADRID OFFICE**

Avenida Ciudad de Barcelona 81,1 Planta 28007, Madrid, SPAIN  
Phone: 01-552-4160 Fax: 01-552-4161

**TOPCON SCANDINAVIA A.B.**

Industrivägen 4 P.O.Box 2140 43302 Sävedalen SWEDEN  
Phone: 031-261250 Fax: 031-268607

**TOPCON SINGAPORE PTE.LTD.**

Alexandra District Park, Block 4, #05-15, Pasir Panjang Road, SINGAPORE 116491  
Phone: 2780222 Fax: 2733540 Telex: 61121 TOPSIN

**TOPCON AUSTRALIA PTY.LTD.**

408 Victoria Road, Gladesville, NSW 2111, AUSTRALIA  
Phone: 02-9617-4656 Fax: 02-9617-4654

**TOPCON INSTRUMENTS (MALAYSIA) SDN. BHD.**

Lot 226 Jalan Negara 2, Pusat Bandar Taman Melawati, Taman Melawati,  
53100 Kuala Lumpur, MALAYSIA  
Phone: 03-4079801 Fax: 03-4079796

**TOPCON INSTRUMENTS (THAILAND) CO.,LTD.**

7th Floor, Thai Virawat Building 86/1, Krungthoburi Road,  
Banglumpoo-Lang, Klongsan Bangkok 10500 THAILAND  
Phone: 662-860-7801-6 Fax: 662-560-7807-8

**TOPCON KOREA CORPORATION**

Hyobong Bldg., 1-1306, Seocho-Dong, Seocho-Gu, Seoul, KOREA  
Phone: 02-557-9231/2 Fax: 02-556-1928

**TOPCON OPTICAL(H.K.) LTD.**

2-4/F Meaco Industrial Bldg., No.53-55 Au Pui Wan Street, Fo Tan Road, Sha Tin,  
N.T. HONG KONG  
Phone: 26901328 Fax: 26910264

**TOPCON CORPORATION BEIJING OFFICE**

Room NO.962 Poly Plaza Building, 14 Dongzhimen  
Nandajie, Dongcheng District, Beijing, 100027 CHINA  
Phone: 10-6501-4191-2 Fax: 10-6501-4190

**TOPCON CORPORATION BEIRUT OFFICE**

P.O. BOX 70-1002 Anielias, BEIRUT-LEBANON  
Phone: 961-1-521119 Fax: 961-1-521119

**TOPCON CORPORATION DUBAI OFFICE**

Office No.102, KHALAF RASHID AL NAYLI BLDG., Deira, Dubai, UAE  
Phone: 971-4-695511 Fax: 971-4-695522